



Pourer „Lotus“
PACKSYS

Pourer „Lotus“



Application advantages:

- Clean and easy pouring
- 100% dropleless dosing
- Single-hand use
- Senior-friendly
- Tamper evidence
- Simple and intuitive handling
- One piece system

Range of application:

Pharmacie:

Cough syrup, syrups, fever reliefs, skin oil (e.g. Neurodermatitis), liquid pain medications, emulsions / suspensions / solutions for ingestion

Cosmetics:

Nurishing oils (body oils, skin oils), bath additives, sun oils, baby oils, low-viscosity lotions

Chemistry and technology:

Detergents and cleansing extract, colours, fertiliser, additives in the water (e.g. for aquaria)

Range of application:

- Pharmacy
- Cosmetics
- Food and beverage
- Nutritional supplements
- Chemistry and technology

Development status:

- Series production from October 2020
- Sampling from March 2020

Equipment:

- Neck finish PP28
- Glass bottle
- Plastic bottle (HDPE, PET, etc.)
- Flip top closure
- Tamper-evident

Characteristics:

- Clean and easy pouring
- Simple and intuitive handling
- One piece system
- Single-hand use

Filling:

- Screwing



pour



drop



One piece
system


Pourer

Tamper-
evident

Benefits for customers and manufacturer:

- Pharmaconform material
- Double tamper-evident
- One-piece system in the colour of your choice
- Customer logo on the hinged cover top side possible
- Fewer filling problems (no missing dosing insert in the closure)
- Large variety of bottles in glass and plastics with PP28 neck finish
- Processing on standard filling plant
- Ideal for pouring different liquids from low to medium viscosity
- At high internal pressure, for example by fermenting liquids, the closure opens before the bursting of the bottle.
- On request combinable with a dosing cup





PACKSYS

PACKSYS GmbH, Alte Weilheimer Str. 2-4, D-82340 Feldafing, Germany
phone +49(0)8157/9251-0, info@packsys.de, www.packsys.de